


# FERBERITE (Iron Wolframite mineral, $\text{FeWO}_4$ )

MILESTONE ultraWAVE 3 METHOD #GEO-03-05



## I SUMMARY

This method provides the acid digestion of the sample in a Single Reactor Chamber (SRC) using temperature control and microwave heating for the metal determination by spectroscopic methods.

## I INSTRUMENTATION FOR DIGESTION



Milestone ultraWAVE 3 with internal temperature and pressure control in all vessels, built-in touchscreen with easyCONTROL software installed.

## I RACK AND VIALS TYPE

| Rack   | Vial |
|--------|------|
| 20 pos | PTFE |

## I SAMPLE AMOUNT AND REAGENTS

| Rack   | Sample amount | Reagents  |
|--------|---------------|---|
| 20 pos | 0,1 g         | 2 mL HCl,<br>2 mL $\text{HNO}_3$ ,<br>2 mL $\text{H}_3\text{PO}_4$ ,<br>1 mL HF |

## I BASE LOAD

120 mL of DI  $\text{H}_2\text{O}$ , 5 mL of  $\text{HNO}_3$  67%

## I TYPE OF GAS AND STARTING PRESSURE

| Gas Type     | Load pressure |
|--------------|---------------|
| Nitrogen gas | 40 bar        |
| Argon Gas    | 40 bar        |

Milestone suggests gas purity of 4.7. For trace metal analysis the purity of 5.0 is highly recommended

## I MICROWAVE PROGRAM

| Step | Time (min) | Power (W) | T1 (°C) | T2 (°C) | P (bar) |
|------|------------|-----------|---------|---------|---------|
| 1    | 00:15:00   | 1500      | 230     | 60      | 100     |
| 2    | 00:15:00   | 1500      | 230     | 60      | 100     |

## I NOTES

- For details on the operating steps with ultraWAVE 3, see the user manual and video tutorials on Milestone Connect.
- Cleaning /maintenance procedure must be executed following the user manual, video tutorial and installer instructions.
- This procedure is only a guideline and it may need to be modified or changed to obtain the required results on your sample.